

Essilor Ideal® Availability

EFFECTIVE DATE: March 2012




Essilor Ideal (17mm MFH)

Materials	Power Range	Cylinder	Add Range	Base Curve
1.50 Clear	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	2.00, 3.50, 4.25, 6.25, 7.25
1.50 Transitions® VI G/B	+6.00 to -8.00	up to -6.00D	0.75 to 3.50	2.25, 4.25, 5.25, 8.25
1.50 Transitions® XTRActive™	+6.00 to -8.00	up to -6.00D	0.75 to 3.50	2.25, 4.25, 5.25, 8.25
Coming Soon 1.50 Transitions® Vantage™	+6.00 to -8.00	up to -6.00D	0.75 to 3.50	2.25, 4.25, 5.25, 8.25
1.50 Xperio® G/B	+6.00 to -10.00	up to -4.00D	0.75 to 3.50	2.00, 3.25, 4.25, 6.25, 8.25
TREXA™ Clear	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	2.00, 3.50, 4.25, 6.25, 7.25
TREXA™ Transitions® VI G/B	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	2.00, 3.50, 4.25, 6.25, 7.25
Airwear® Clear	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Transitions® VI G/B	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Transitions® XTRActive™	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Transitions® Vantage™	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Xperio® G/B	+8.00 to -10.00	up to -6.00D	0.75 to 3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50
Thin&Lite® 1.60 Clear	+8.00 to -12.00	up to -6.00D	0.75 to 3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50
Thin&Lite® 1.60 Transitions® G/B	+8.00 to -12.00	up to -6.00D	0.75 to 3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50
Thin&Lite 1.67 Clear	+9.00 to -12.00	up to -6.00D	0.75 to 3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00
Thin&Lite 1.67 Transitions G/B	+9.00 to -12.00	up to -6.00D	0.75 to 3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00
Thin&Lite® 1.67 Xperio® G/B	+9.00 to -12.00	up to -4.00D	0.75 to 3.50	1.75, 2.75, 3.75, 5.25, 6.75, 8.00

Essilor Ideal Short™ (14mm MFH)

Materials	Power Range	Cylinder	Add Range	Base Curve
1.50 Clear	+6.00 to -10.00	up to -6.00D	0.75 to 3.00	2.00, 3.50, 4.25, 6.25, 7.25
1.50 Transitions® VI G/B	+6.00 to -8.00	up to -6.00D	0.75 to 3.00	2.25, 4.25, 5.25, 8.25
1.50 Transitions® XTRActive™	+6.00 to -8.00	up to -6.00D	0.75 to 3.00	2.25, 4.25, 5.25, 8.25
Coming Soon 1.50 Transitions® Vantage™	+6.00 to -8.00	up to -6.00D	0.75 to 3.00	2.25, 4.25, 5.25, 8.25
1.50 Xperio® G/B	+6.00 to -10.00	up to -4.00D	0.75 to 3.00	2.00, 3.25, 4.25, 6.25, 8.25
TREXA™ Clear	+6.00 to -10.00	up to -6.00D	0.75 to 3.00	2.00, 3.50, 4.25, 6.25, 7.25
TREXA™ Transitions® VI G/B	+6.00 to -10.00	up to -6.00D	0.75 to 3.00	2.00, 3.50, 4.25, 6.25, 7.25
Airwear® Clear	+6.00 to -10.00	up to -6.00D	0.75 to 3.00	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Transitions® VI G/B	+6.00 to -10.00	up to -6.00D	0.75 to 3.00	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Transitions® XTRActive™	+6.00 to -10.00	up to -6.00D	0.75 to 3.00	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Transitions® Vantage™	+6.00 to -10.00	up to -6.00D	0.75 to 3.00	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Xperio® G/B	+8.00 to -10.00	up to -6.00D	0.75 to 3.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50
Thin&Lite® 1.60 Clear	+8.00 to -12.00	up to -6.00D	0.75 to 3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50
Thin&Lite® 1.60 Transitions® G/B	+8.00 to -12.00	up to -6.00D	0.75 to 3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50
Thin&Lite 1.67 Clear	+9.00 to -12.00	up to -6.00D	0.75 to 3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00
Thin&Lite 1.67 Transitions G/B	+9.00 to -12.00	up to -6.00D	0.75 to 3.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00
Thin&Lite® 1.67 Xperio® G/B	+9.00 to -12.00	up to -4.00D	0.75 to 3.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00

G/B=Available in both Gray and Brown



Essilor Ideal® Availability

EFFECTIVE DATE: March 2012

Essilor Ideal Advanced™ (17mm MFH)

Materials	Power Range	Cylinder	Add Range	Base Curve
1.50 Clear	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	2.00, 3.50, 4.25, 6.25, 7.25
1.50 Transitions® VI G/B	+6.00 to -8.00	up to -6.00D	0.75 to 3.50	2.25, 4.25, 5.25, 8.25
1.50 Transitions® XTRActive™	+6.00 to -8.00	up to -6.00D	0.75 to 3.50	2.25, 4.25, 5.25, 8.25
Coming Soon 1.50 Transitions® Vantage™	+6.00 to -8.00	up to -6.00D	0.75 to 3.50	2.25, 4.25, 5.25, 8.25
1.50 Xperio® G/B	+6.00 to -10.00	up to -4.00D	0.75 to 3.50	2.00, 3.25, 4.25, 6.25, 8.25
TREXA™ Clear	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	2.00, 3.50, 4.25, 6.25, 7.25
TREXA™ Transitions® VI G/B	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	2.00, 3.50, 4.25, 6.25, 7.25
Airwear® Clear	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Transitions® VI G/B	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Transitions® XTRActive™	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Transitions® Vantage™	+6.00 to -10.00	up to -6.00D	0.75 to 3.50	1.50, 3.25, 4.25, 5.50, 7.00
Airwear® Xperio® G/B	+8.00 to -10.00	up to -6.00D	0.75 to 3.50	1.50, 3.25, 4.25, 5.50, 7.00, 8.50
Thin&Lite® 1.60 Clear	+8.00 to -12.00	up to -6.00D	0.75 to 3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50
Thin&Lite® 1.60 Transitions® G/B	+8.00 to -12.00	up to -6.00D	0.75 to 3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.50
Thin&Lite 1.67 Clear	+9.00 to -12.00	up to -6.00D	0.75 to 3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00
Thin&Lite 1.67 Transitions G/B	+9.00 to -12.00	up to -6.00D	0.75 to 3.50	1.75, 2.75, 3.75, 6.00, 6.75, 8.00
Thin&Lite® 1.67 Xperio® G/B	+9.00 to -12.00	up to -4.00D	0.75 to 3.50	1.75, 2.75, 3.75, 5.25, 6.75, 8.00

Essilor Ideal Advanced Wrap™

Materials	Power Range	Cylinder	Add Range	Base Curve
1.50 Clear	-3.50 to +3.00	up to -3.00D	0.75 to 3.00	6.25, 8.00
1.50 Transitions VI G/B	-4.00 to +3.00	up to -4.00D	0.75 to 3.50	6.25, 8.25
1.50 Xperio G/B	-4.00 to +3.00	up to -3.50D	0.75 to 3.00	6.25, 8.25
Airwear Clear	-4.00 to +3.00	up to -4.00D	0.75 to 3.50	6.50 (HC only), 8.50
Airwear Transitions VI G/B	-4.00 to +3.00	up to -4.00D	0.75 to 3.50	6.50, 8.50
Airwear Xperio G/B	-4.00 to +3.00	up to -4.00D	0.75 to 3.50	6.50, 8.50
Thin&Lite 1.67 Clear	-5.00 to +3.00	up to -4.00D	0.75 to 3.50	6.00, 8.00 (UC only)
Thin&Lite 1.67 Transitions VI G/B	-5.00 to +3.00	up to -4.00D	0.75 to 3.50	6.00, 8.00

G/B=Available in both Gray and Brown